


# 48-FPC-2-2 dc-xxx Mod01

Fiber Port Cluster 2 → 2 dichroic, with shutters



#### **FEATURES**

Fiber Port Cluster for one input source

- Configuration 2 → 2 dc
- Electro-magnetic shutters at the two input ports
- Electro-magnetic shutters at the two output ports
- Highly efficient coupling into polarizationmaintaining fiber cables
- Adjustable splitting ratio
- Compact, rugged, transportable and sealed optomechanical units
- Fully fiber-coupled
- Very high long-term stability, efficiency and reproducability

## **DESCRIPTION**

This Fiber Port Clusters  $2 \rightarrow 2$  dc Mod01 is a compact opto-mechanical unit that combines two fiber-coupled sources with differing wavelengths and then splits the combined radiation into 2 output fiber cables with high efficiency and variable splitting ratio. Electro-magnetic shutter is placed at each of the two input ports and output ports.

#### **Optical Setup**

The input port is fiber-coupled to a <u>PM fiber cable</u>. A polarizer defines the input polarization which is necessary for a long term stable splitting ratio.

A photo diode right after the input port allows for a continuous monitoring of the radiation.

Subsequently, the radiation splitting is achieved by using a rotary half-wave plate in combination with a polarization beam splitter. By use of the rotary half-wave plate, almost any desired splitting ratio can be achieved.

At the output ports further polarizers are placed in order to define the polarization at output of the system.

#### **Fiber Couplers**

A fundamental component of a Fiber Port Cluster is the <u>Laser Beam Coupler</u>, which is the input into the opto-mechanical unit collimating the input radiation and, finally, couples the radiation back into the polarization-maintaining fiber cables. The stability of the total Fiber Port Cluster is determined by the <u>stability</u> of the laser beam coupler.




#### **Electro-magnetic Shutters**

An electro-magnetic shutter is placed at each of the two input ports and output ports. For a description and further technical specifications, please refer to the electro-magnetic shutter type <u>48EMS-6</u>.

#### How to order

For a detailed quotation please additionally specify

- Wavelengths In 1 and In 2
- Cable lengths
- Connector types



# **TECHNICAL DATA**

48-FPC-2-2\_dc-xxx\_Mod01

Order code	48-FPC-2-2_dc-xxx_Mod01
Configuration	2 → 2 dc
	electro-magnetic shutters
Wavelengths*	767 + 780 nm
Fiber type	polarization-maintaining
Connector type	FC APC (standard)
Cable lengths	customer-specific
Wave plate type	low-order
Power monitor	BPX-61 (SMA)
Transmission	≥ 70 % @ 767 nm
	≥ 60 % @ 780 nm
Polarization Extinction Ratio	≥ 23 dB @ 767 + 780 nm
Balancing	better 3 %



\* Different wavelength combinations on request

# **TECHNOTES**

- Article Fiber Port Cluster Rugged, modular and fiber coupled beam splitting and combining units
- Connecting multicube assemblies to a base plate How to connect the self-supporting multicube system

### **DOWNLOADS**



980129090615.pdf (Dimensional drawing)



Article Cluster.pdf (Technote)

### **RELATED PRODUCTS**

**FIBER COLLIMATOR** 

60FC-Q

Fiber Collimator for collimating large beam diameters

and with integrated quarter-wave plate

**POLARIZATION** 

**ANALYZER SK010PA** 

Measurement tool for coupling into polarization-

maintaining fiber cables

**FIBER COLLIMATOR** 

**SERIES 60FC-SF** 

Fiber Collimator/Fiber Coupler with super-fine thread

**ELECTRO-MAGNETIC SHUTTER 48EMS-6** 



This is a printout of the page https://sukhamburg.com/products/details/48-FPC-2-2 dc-xxx Mod01 from 4/29/2024

# **CONTACT**

For more information please contact: Schäfter + Kirchhoff GmbH Kieler Str. 212 22525 Hamburg Germany

Tel: +49 40 85 39 97-0 Fax: +49 40 85 39 97-79

info@sukhamburg.de www.sukhamburg.com

# **LEGAL NOTICE**

Copyright 2020 Schäfter+Kirchhoff GmbH. All rights reserved.

Text, image, graphic, sound, video and animation files and their arrangement on Schäfter+Kirchhoff GmbH webpages are protected by copyright and other protective laws. The content may not be copied for commercial use or reproduced, modified or used on other websites. [more]